

Immunofluorescence test (tube method)

1. For each sample to be tested (no. of tests plus controls x no. of donors) aliquot 50uL platelet suspension into 'FACS tubes'(Falcon 2052)
2. Add 50uL test serum/plasma. Include 2 negative controls and 1 weak positive control which has been calibrated against the National anti-HPA-1a standard (NIBSC code 93/710)
3. Incubate at 37°C for 30 min.
4. Add 4ml PBS/EDTA buffer and spin for 5 min at 1,500g in benchtop centrifuge.
5. Decant supernatant, re-suspend cells, add 4mL PBS/EDTA wash buffer and spin for 5 min at 1,500g.
6. repeat step 5.
7. Decant supernatant, re-suspend cells and add 50uL FITC-conjugated anti-human Ig (G, M or G+M+A, dilution determined by previous experiments). Make sure cells are fully re-suspended.
8. Incubate at 22°C (room temperature) for 30 mins in the dark.
9. Repeat step 4. and decant supernatant.
10. Re-suspend cells in 1mL PBS/EDTA buffer.
11. Store in the dark until FACS analysis.

Immunofluorescence test (microplate method)

1. For each sample to be tested (no. of tests plus controls x no. of donors) aliquot 50uL platelet suspension into U-well microplate and spin at 900g for 3 min.
2. Discard supernatant and re-suspend cells in 50uL test serum/plasma.
3. Include 2 negative controls and 1 weak positive control which has been calibrated against the National anti-HPA-1a standard (NIBSC code 93/710)
4. Incubate at 37°C for 30 min.
5. Add 200ul PBS/EDTA buffer and spin for 3 min at 900g in microplate centrifuge.
6. Decant supernatant, re-suspend cells using microplate shaker at maximum speed for 20 secs, add 200uL PBS/EDTA buffer and spin for 3 min at 900g.
7. repeat step 5.

8. Decant supernatant, re-suspend cells using microplate shaker and add 50uL FITC-conjugated anti-human Ig (G, M or G+M+A, dilution determined by previous experiments). Make sure cells are fully re-suspended using microplate shaker or pipette.
9. Incubate at 22°C (room temperature) for 30 mins in the dark.
10. Repeat step 5. and decant supernatant.
11. Re-suspend cells in 100uL PBS/EDTA wash buffer and transfer to numbered 'FACS tubes' containing 1mL PBS/EDTA wash buffer.
12. Store in the dark until FACS analysis.

N.B If using solid phase technique; use 10uL volumes of platelet suspension, test sera, and FITC-anti-IgG and wash 10 times at steps 6. and 10. If using visual microscopy to read tests carried out in microplates or tubes then cells should be washed twice more after step 10.