**CE Marked Material**  
**QCRHBcIgMQC1 - IgM Anti-HBc Quality Control Reagent Sample**  
**1**

**NIBSC code: 07/B498-xxx**  
**Instructions for use**  
**(Version 9.0, Dated 29/10/2019)**

**This material is an 'Annex II List A' IVD and complies with the requirements of the "EU in vitro diagnostic medical device directive 98/79/EC".**

### 1. INTENDED USE

**This product is CE marked for use as an IVD within the EU member states and EEA countries. In all other territories this product can be used for research purposes only.**

IgM anti-HBc QC1 is intended for use in the internal laboratory quality control of immunoassays that detect IgM antibodies to hepatitis B core antigen. The IgM anti-HBc QC1 should be included in each run as part of a continuing quality control programme to monitor the performance of the assay. Data obtained with the IgM anti-HBc QC1 can be used to construct quality control charts that can be visually monitored each time the assay is run, to check for consistency of performance of the assay. Examples of how these charts are constructed and used have been described elsewhere<sup>1</sup>. IgM anti-HBc QC1 is NOT INTENDED TO BE USED TO COMPARE THE SENSITIVITY OF PARTICULAR ASSAYS.

### 2. CAUTION

**This preparation is not for administration to humans or animals in the human food chain.**

The IgM anti-HBc QC1 has been prepared from a pool of IgM anti-HBc reactive defibrinated donations, repeatedly reactive in commercial EIA kits. The reactive donation used to prepare IgM anti-HBc QC1 was non-reactive for HBsAg, anti-HIV and anti-HCV using commercial EIA kits. The reactive serum was then diluted in a pool of defibrinated human plasma samples non-reactive for IgM anti-HBc. These samples were also non-reactive for HBsAg, anti-HCV and anti-HIV using commercial EIA kits. Bronidox® was added to a concentration of 0.05%(w/v) as a preservative. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

### 3. UNITAGE

**Table 1** gives a summary of the results obtained for IgM anti-HBc QC1 **07/B498** These results are intended only as a guide to the approximate levels of reactivity to be expected, and may not be exactly reproduced in other laboratories. In each case, at a minimum, three samples of IgM anti-HBc QC1 were tested on two separate occasions. The results are expressed as the ratio of mean optical density or other measurement of the IgM anti-HBc response of the QC1 sample, to the kit manufacturer's calculated cut-off.

### 4. CONTENTS

Country of origin of biological material: United Kingdom.  
Ready-to-use reagent  
REF QCRHBcIgMQC1 1x1mL Sarsdedt Bottles  
Defibrinated Plasma 1mL  
Bronidox® (Sigma-Aldrich) 0.05% (w/v)

### 5. STORAGE

- Reagents are to be kept at 2-8°C upon receipt
- Reagents may be stored at 2-8°C until use by date
- Reagents should be divided into measured sub-aliquots and stored below -20°C to avoid freeze/thaw cycles.

National Institute for Biological Standards and Control,  
Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, [nibsc.org](http://nibsc.org)  
WHO International Laboratory for Biological Standards,  
UK Official Medicines Control Laboratory

- When thawed for use, store at 2-8°C. Once thawed, use within one month and do not refreeze
- Ensure all containers are properly sealed to avoid drying out of the reagent
- Avoid microbial contamination of this product as this may alter product performance
- Avoid excessively high temperatures or humidity

### 6. DIRECTIONS FOR OPENING

Vials have a screw cap; an internal stopper may also be present. The cap should be removed by turning anti-clockwise. Care should be taken to prevent loss of the contents. Please note: If a stopper is present on removal of the cap, the stopper should remain in the vial or be removed with the cap.

### 7. USE OF MATERIAL

1. Use of this reagent is to be restricted to trained laboratory staff only
2. Use suitable (latex/nitrile) gloves and eye/skin protection
3. Include reagent as a normal sample in routine work list
4. Allow reagent to reach room temperature before use
5. Plot reagent result on a QC chart to monitor performance

### 8. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

Stability study analysis has been performed on this reagent to assess stability. The expiry date is shown on the label.

### 9. REFERENCES

1. Levey, S. and Jennings, E.R. (1950) The use of control charts in clinical laboratories. *Am.J.Clin.Pathol.* 20, 1059-1066

### 10. ACKNOWLEDGEMENTS

### 11. FURTHER INFORMATION

Further information can be obtained as follows:  
This material: [enquiries@nibsc.org](mailto:enquiries@nibsc.org)  
WHO Biological Standards:  
<http://www.who.int/biologicals/en/>  
JCTLM Higher order reference materials:  
<http://www.bipm.org/en/committees/jc/jctlm/>  
Derivation of International Units:  
[http://www.nibsc.org/standardisation/international\\_standards.aspx](http://www.nibsc.org/standardisation/international_standards.aspx)  
Ordering standards from NIBSC:  
<http://www.nibsc.org/products/ordering.aspx>  
NIBSC Terms & Conditions:  
[http://www.nibsc.org/terms\\_and\\_conditions.aspx](http://www.nibsc.org/terms_and_conditions.aspx)

### 12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

### 13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

### 14. MATERIAL SAFETY SHEET

Physical and Chemical properties	
<b>Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified</b>	Corrosive: No
Physical appearance: Liquid	
Stable: Yes	Oxidising: No
Hygroscopic: No	Irritant: No
Flammable: No	Handling: See caution, Section 2
Other (specify):	
Toxicological properties	
Effects of inhalation:	Not established, avoid inhalation
Effects of ingestion:	Not established, avoid ingestion
Effects of skin absorption:	Not established, avoid contact with skin
Suggested First Aid	
Inhalation:	Seek medical advice
Ingestion:	Seek medical advice
Contact with eyes:	Wash with copious amounts of water. Seek medical advice
Contact with skin:	Wash thoroughly with water.
Action on Spillage and Method of Disposal	
Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.	

#### 15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at [http://www.nibsc.org/About\\_Us/Terms\\_and\\_Conditions.aspx](http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx) or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

#### 16. INFORMATION FOR CUSTOMS USE ONLY

<b>Country of origin for customs purposes*:</b> United Kingdom * Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
<b>Net weight:</b> 1.5 g
<b>Toxicity Statement:</b> Toxicity not assessed
<b>Veterinary certificate or other statement</b> if applicable. <b>Attached:</b> No

**TABLE 1:** Results obtained for **IgM Anti-HBc QC1** (Lot Number **07/B498**) using the following EIA kits.

EIA KIT	Method Options	Test to Cut-off Ratio	
		Mean	SD (n-1)
<b>Monolisa Anti-HBc IgM</b> Manufacturer: Bio-Rad Catalogue number: 72382 Lot number: 7A0021	Standard Protocol	2.07	0.08
<b>AxSYM Core-M *</b> Manufacturer: Abbott Diagnostics Catalogue number: 7A44-20 Lot number: 48477LUOO & 47166LUOO	Automated	1.24	0.04
<b>Bioelisa HBc IgM</b> Manufacturer: Biokit Catalogue number: 3000-1103 Lot number: K-3506	Standard Protocol	2.90	0.24
<b>Murex Anti-HBc IgM</b> Manufacturer: Abbott Diagnostics Catalogue number: 9E25-01 Lot number: J267510	Standard Protocol - Qualitative	1.52	0.10
<b>Bayer Centaur ~</b> Manufacturer: ADVIA Catalogue number: 00504619 Lot number: 9382213	Automated	1.61	0.08

\* Tests performed at Northwick Park Hospital and The William Harvey Hospital

~ Tests performed at St. Mary's Hospital

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<http://www.bipm.org/en/committees/jc/jctlm/>

Derivazione delle Unità Internazionali:  
[http://www.who.int/biologicals/reference\\_preparations/en/](http://www.who.int/biologicals/reference_preparations/en/)  
 Ordine standards da NIBSC:  
[http://www.nibsc.ac.uk/products/ordering\\_information/frequently\\_asked\\_questions.aspx](http://www.nibsc.ac.uk/products/ordering_information/frequently_asked_questions.aspx)  
 NIBSC Termini e condizioni:  
[http://www.nibsc.ac.uk/terms\\_and\\_conditions.aspx](http://www.nibsc.ac.uk/terms_and_conditions.aspx)

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