CE Marked Material
Anti Rubella Serum, Human
NIBSC code: 67/182
Instructions for use
(Version 10.0, Dated 31/03/2020)

This material is an 'Annex II List B' IVD and complies with the requirements of the "EU in vitro diagnostic medical device directive 98/79/EC".

1. INTENDED USE

This product is CE marked for use as an IVD within the UK, EU member states and EEA countries. In all other territories this product can be used for research purposes only.

The product is intended for use in the field of in vitro diagnostics, in the detection and quantitation of rubella infection.

The Standard is not intended for administration to humans but is supplied to 'professional' users only.

The 2nd British Standard for Anti-Rubella Serum, Human (67/182) was established by the National Biological Standards Board (NBSB) in 1986.

It was initially and indirectly quantitated in 1974, being one of the preparations (Serum B 67/182) used in an international collaborative study in the calibration of the 1st British Standard for Anti-Rubella Serum, Human (69/60) against the 2nd WHO International Reference Preparation of Anti-Rubella Serum, Human¹.

In 1984 when stocks of the 1st British Standard were running low an obvious replacement candidate was the human anti-rubella serum 67/182. It had a significant advantage since it has been calibrated in the 1974 collaborative study; it had given every indication of being a stable preparation over the intervening time and finally it was available in sufficient numbers. This preparation of human anti-rubella serum 67/182 was then sent out in a collaborative study to five laboratories²

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV1 and 2 and HCV RNA³. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

Each ampoule contains the following: **80 International Units (IU)**

4. CONTENTS

Country of origin of biological material: United Kingdom.

The 2nd British Standard for Anti-Rubella Serum, 67/182 consists of ampoules each containing equal quantities of the residue after freeze-drying of approximately 1ml of pooled hyperimmune human sera.


5. STORAGE

This product should be stored -20°C or below.

6. DIRECTIONS FOR OPENING

Tap the ampoule gently to collect the material at the bottom (labelled) end. Ensure ampoule is scored all round at the narrow part of the neck, with a diamond or tungsten carbide tipped glass knife file or other suitable implement before attempting to open. Place the ampoule in the ampoule opener, positioning the score at position 'A'; shown in the diagram below. Surround the ampoule with cloth or

layers of tissue paper. Grip the ampoule and holder in the hand and squeeze at point 'B'. The ampoule will snap open. Take care to avoid cuts and projectile glass fragments that enter eyes. Take care that no material is lost from the ampoule and that no glass falls into the ampoule.

Side view of ampoule opening device containing an ampoule positioned ready to open. 'A' is the score mark and 'B' the point of applied pressure.

7. USE OF MATERIAL

No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution

Dissolve the total contents of the ampoule by washing out both parts with distilled water or a solution of a suitable buffer and adjust the solution so obtained to a known volume. The solution will contain 80 INTERNATIONAL UNITS. The contents of the ampoule are extremely hygroscopic and no attempt should be made to weigh out any of the freeze dried material.

8. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

Once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use.

In 1971 a number of unopened ampoules of the human anti-rubella serum (67/182) were placed at elevated temperatures in an accelerated degradation study. In 1985, ampoules which had been kept at +4°C for 14 years were assayed against other ampoules of the bulk of the batch which had been stored at -20°C (the normal holding temperature for standards) for the same time span, and a small but statistically significant (P<0.05) loss of potency was detected in the ampoules held at +4°C.

In 2005 a further ampoules which had been kept at +4°C were assayed against ampoules stored at -20°C. These ampoules had been stored at elevated temperatures for 34 years. A slight loss of potency relative to the -20°C samples was observed but the loss was not as great as for the assays from 1985.

This suggests that after an initial slight drop, the +4°C samples have remained stable for 20 years with no detectable loss. In 2015 a further study demonstrated that the material has exhibited no loss in potency of the +4°C samples stored for 45 years. On the basis of this stability data the expiry date assigned to the product is 2025-09-30.

9. REFERENCES

1. Clarke, M. et al (1975), Journal of Biological Standardization, 3, 151.
2. Unpublished document (1986), a copy of which may be obtained by making application to this Institute.

10. ACKNOWLEDGEMENTS

The bulk material consisted of pooled hyperimmune serum obtained from blood donors and generously provided by the National Blood Transfusion Service (NBS).

EC REP: Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta.

11. FURTHER INFORMATION

Further information can be obtained as follows:

This material: enquiries@nibsc.org
 WHO Biological Standards:
<http://www.who.int/biologicals/en/>
 JCTLM Higher order reference materials:
<http://www.bipm.org/en/committees/jc/jctlm/>
 Derivation of International Units:
http://www.nibsc.org/standardisation/international_standards.aspx
 Ordering standards from NIBSC:
<http://www.nibsc.org/products/ordering.aspx>
 NIBSC Terms & Conditions:
http://www.nibsc.org/terms_and_conditions.aspx

reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom
* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
Net weight: 0.1g
Toxicity Statement: Non-toxic
Veterinary certificate or other statement if applicable.
Attached: No

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties	
Physical appearance: Freeze-dried powder	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: Yes	Irritant: No
Flammable: No	Handling: See caution, Section 2
Other (specify):	Contains material of human origin
Toxicological properties	
Effects of inhalation:	Not established, avoid inhalation
Effects of ingestion:	Not established, avoid ingestion
Effects of skin absorption:	Not established, avoid contact with skin
Suggested First Aid	
Inhalation:	Seek medical advice
Ingestion:	Seek medical advice
Contact with eyes:	Wash with copious amounts of water. Seek medical advice
Contact with skin:	Wash thoroughly with water.
Action on Spillage and Method of Disposal	
Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.	

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by

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