

Data Sheet

| | |
|------------------------------|---|
| PRODUCT NAME: | PepMix ESAT-6 |
| PRODUCT CODE: | PM-MYCTU-ESAT6 |
| REPOSITORY REFERENCE: | 100034 |
| DESCRIPTION: | 21 peptides, 15meric; approx. 25 µg / peptide |
| SOURCE: | esxA 6 kDa early secretory antigenic target ESAT-6 Uniprot ID : P9WNK7 |
| PRESENTATION: | The peptides of this product are supplied as trifluoroacetate salts. |
| PURITY (HPLC): | >70% (determined at 220 nm; individual peptides) |
| STORAGE: | Store at - 20°C |
| SOLUBILITY | |
| RECOMMENDATION: | Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system. |
| CONTRIBUTOR: | JPT Peptide Technologies GmbH, Volmerstrasse 5, 12489 Berlin, Germany |

ACKNOWLEDGEMENTS

Publications should acknowledge the contributor and the Centre for AIDS Reagents. Acknowledgments should read: "The *Name of Reagent (Repository Number)* was obtained from the Centre for AIDS Reagents, NIBSC, UK, supported by EURIPRED (EC FP7 INFRASTRUCTURES-2012 - INFRA-2012-1.1.5.: Grant Number 31266). www.euripred.eu/